

Single-cell RNA-sequencing (scRNA-seq): knowing the in and outs of the data generated

Ken Lau, Assistant Professor of Cell and Developmental Biology

ken.s.lau@vanderbilt.edu

CQS Summer Academy (8/13/2018)

<http://www.mc.vanderbilt.edu/vumcdept/cellbio/laulab/index.html>

Twitter: @KenLauLab



About me and single cells

- Started lab at Vanderbilt in 2013 with focus on single-cell biology of the gut (IBD and colon cancer)
- Multiplex imaging, CyTOF, scRNA-seq
- Training at Toronto/MIT/Harvard on multivariate analysis, mathematical modeling, and tissue systems
- inDrop in lab since August 2016 (first 1cell customer outside of Boston), > 50 samples ran so far > 500 000 cells sequenced; we have two systems



Outline

- Introduction to scRNA-seq techniques
- Discussion on scRNA-seq data issues
- Brief Python introduction

https://github.com/KenLauLab/Discovery_Oriented_Data_Science

scRNA-seq techniques

Bulk data versus single-cell data






Bulk

VS.




Single-cell

Bulk RNA-seq



1 sample x 25,000 Genes

Single-cell RNA-seq



scRNAseq protocols

Cell Encapsulation techniques

- Droplet-Based
- Well-Based
- Microfluidic capture (Fluidigm C1)

Lysis and RT

- Coupled – requires balanced mix
- Uncoupled – enables more aggressive lysis

RNA capture strategies

- Poly dT priming
- Targeting / enrichment


Indexing strategies

- During capture/RT (typically per **cell** indexes - barcodes)
- After RT (typically per **well** indexes)


Amplification strategies

- PCR vs IVT



Typical logic for a scRNA-seq experiment




Single cell suspension***



Single-cell encapsulation/
Library preparation



Sequencing and alignment
(Bioinformatics I)




Well-known methods to isolate single cells

Laser-capture microdissection


(i) Cell selection



(ii) Laser sectioning



(iii) Cell transfer on a membrane



Fluorescence-activated cell sorting






Plate-Based scRNAseq





- Isolate RNA, label transcripts using barcoded RT primers (3' seq) or through template switching library prep (enables full length)





Relative to other platforms:

- $\sim 10\mu\text{l}$ / cell, < 1000 cells (higher volume, lower throughput)
- Deeper sequencing possible - flexibility
- Number of wells become limiting (doublet rate vs. cost)

Microfluidic capture scRNAseq (Fluidigm C1)



Capture chips are built to accommodate specific ranges of cell size




Cells

Stanford course

Bead-based capture: Immobilized sets of indexed primers

- Each bead is coated with primers containing a barcode unique to that bead – (index for each cell)



Barcode



UMI

TTTT

- And a Unique Molecular Identifier (UMI) that uniquely tags each primer – (index for each transcript)

Droplet-based encapsulation

- Co-encapsulating cells and beads in thousands of 1-5nL droplets
- Beads carry barcoded poly-T primers to capture RNA
- Encapsulation rate follows Poisson distribution
- Excess of “vessels” to minimize doublets



Cells entering/ sec x


Beads/droplet

=

cells captured/ sec

Beads


- inDrop
(1cellBio)



- Drop-seq
(Chemgenes)




- Chromium
(10x Genomics)



Microfluidic droplet encapsulation chip

Beads




Cells


Droplet
formation

Emulsion output


Seq-Well - microwell sequencing (Shalek lab)




Drop-Seq




Seq-Well



Seq-Well: Principle




Nanowell Array



Size Exclusion → ≤ 1 bead per well

Sealing → ~~Cross Contamination~~

Deconvolving the data



Reads with same barcodes collapse into cells

Read with same UMIs collapse into transcript counts

	Cell:	1	2	...	N
Count unique UMIs for each gene in each cell	GENE 1	1	2		14
	GENE 2	4	27		8
	GENE 3	0	0		1
	⋮	⋮	⋮		⋮
	⋮	⋮	⋮		⋮
Create digital expression matrix	GENE M	6	2		0

scRNA-seq Data Exploration (and problems)





Table of genes and barcodes

Genes


[illegible]

Barcodes

Inflection point method for identifying barcodes with real cells



Doublet rate determination




Procedure for the isolation of high viability single-cells from tissues

Keep all reagents on ice and perform all procedures at 4°C, avoid working with overconcentrated tissue/cell solutions


Tissue

*Keep cold and minimize ischemic time

*Mince if needed



Normal colonic epithelium can be isolated through chelation, while other tissues may require direct mechanical processing to achieve 50-500µm fragments. Care should be taken to remove dead cells during washes to maintain viability



The process of isolating tissue fragments should be optimized to the needs of the target cells. This may mean filtering or taking other measures to enrich.

Wash in PBS
Decant small debris

Chelation

10ml DPBS (-Mg/Ca)
3mM EDTA
0.5mM DTT
10mM NAC



~30-60min
Change buffer every 10-15 min


Use forceps to move tissue

Move back to fresh chelation buffer

Shaking

10ml DPBS (-Mg/Ca) in a new tube

2-3 Shakes/sec
~6 inches motion




Wash & inspect crypts

Medimachine

In cold room

Pre-wet 50µm medicon before use
Wash quickly and thoroughly after




*GentleMACS is a potentially less damaging approach we have not yet tested for breaking apart tissue and further dissociating cells.

Wash in PBS
Decant small debris

Washing


DPBS (-Mg/Ca) spin @~300xg for 2.5min



Good pelleting:
(fat or pancreas may still float in first wash)



Sticking to tube;
may need more chelation (EDTA)*



Floating/clumping
typically caused by excess mucin or inviable cells**

Cold Protease Dissociation


2-5ml DPBS (-Mg/Ca) + 2.5mg/ml DNase + 5mg/ml Subtilisin

*Use 2ml for every 50-100µl of pelleted tissue, don't overconcentrate

*Add buffer to frozen subtilisin aliquot to thaw immediately before use

*Pipette with a 1000µl tip every 5-10 min and check for singlets

4-6°C with
gentle motion



*Pellet without EDTA prior to re-suspending in protease

**Optimize spin conditions so that single cells are decanted


**Don't leave tissue in pellet too long, split if too dense

**Re-suspend well, but do not over pipette prior to dissociation


Single-cell isolation for scRNA-seq to minimize dissociation artifacts

- cold protease from *Bacillus licheniformis*, soil bacteria from Himalayan glaciers
- enables tissue preparation on ice (at 4 degrees)

37 degree DNase/Collagenase




Cold Protease (on ice)




- caveat – the efficacy for retrieving all cell types from all tissues unknown

Stressed/dying cells


PCA of mitochondrial gene expression



Contaminant from free floating RNA and leaky cell corpses




Generous threshold



Strict threshold


Amplification bias in scRNA-seq




Non-linear amplification

- Highly expressed transcripts are inappropriately represented and replicated
- Sampling of RNA in a cell 1-10%
- Zero inflated data
- Count data – negative binomial distribution

In vitro transcription results in linear amplification of RNA




Negative binomial versus Gaussian distributions



Gaussian

- Continuous variable
- Symmetric



Negative Binomial

- Discrete variable (counts)
- Asymmetric at small means

Zero inflation in scRNA-seq

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
2	0	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
3	0	0	0	0	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
4	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0
5	0	0	0	0	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
6	0	0	0	0	0	0	0	1	0	0	1	0	0	0	0	0	1	0	0	0	0	0	0	0	0
7	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0
8	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
9	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
10	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
11	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
12	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0
13	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1
14	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
15	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1
16	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0
17	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
18	0	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
19	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
20	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
21	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
22	0	0	0	0	0	0	0	1	0	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
23	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
24	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
25	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
26	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	2	0	0	0	0
27	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
28	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
29	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
30	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
31	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
32	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1
33	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0
34	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
35	0	0	0	0	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
36	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
37	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0

- Only 5% of this table is not 0!

What does this mean?

- Statistical tests to determine differentially expressed genes with a zero-inflated negative binomial (ZINB)
- False negatives (genes that are supposed to be expressed by appear as 0)
- Many genes (columns) have low counts due to shallow sampling of transcripts – this means the data are noisy
- Unreliable variables that need to be processed/filtered out prior to downstream analysis

Feature Selection

Variance selection


- Easiest – select genes that are the most variable

- Variance = genes that are most different across all cells

$$\sigma^2 = \sum \frac{(X - \mu)^2}{N}$$

- Rank genes by top 500 most variable, for example, discard the rest


Highly variable axes = Principal Components



PCA (Principal Component Analysis)

- Principle of PCA is to maximize the Variance of X with the least amount of principal components (latent variables)
- What is variance? Spread of the data, information content, change etc.
- Variance is the covariance of a dataset with itself, i.e. $\text{Var}(X) = \text{Cov}(X,X) \rightarrow \text{Maximize}$
- What are principal components? Linear combinations of original variables – linear transformation


Vectors and projections



Basis Set

$(1,0)$ \longrightarrow
 $(0,1)$ \uparrow


\longrightarrow
Projection




Basis Set

$(7,4)$
 $(-3.5,7)$

PCA as a dimension reduction tool



Selection of principal components



PRINCIPAL COMPONENTS CAN BE ORDERED BY EIGENVALUES (VARIANCE CAPTURED)

Simple summary of this simple feature selection procedure

- Keep most variable genes (over all cells) for downstream analysis, discard rest as noise
- PCA identifies super axes (Principal Components) that are combinations of the original variables that captures the most variance in the data (by eigenvalues)
- Orthogonal – no duplicate or redundant axes – so will only have a few of them